

# CryoFit<sup>®</sup> DMSO

**DMSO for cryopreservation  
(10 ml vial)**



***“CE Certified”***

Cryopreservation of bone marrow, hemopoietic stem cells, human umbilical cord blood and biological cells.

Appropriate for all Healthcare applications!

## SPECIFICATIONS

<b>Sterile</b>	<b>Dimethyl Sulfoxide USP, Ph.Eur.</b>
<b>Endotoxin free</b>	<b>Formula</b> $C_2H_6OS$
<b>Pyrogen free</b>	<b>CAS No</b> 67-68-5
<b>Hygroscopic</b>	<b>EINECS No</b> 200-664-3
<b>Store between 20-30 °C</b>	<b>MW</b> 78.13 g/mol
<b>Protect from strong light</b>	<b>Melting point</b> 18°C
<b>12x10 ml vials</b>	<b>Purity</b> > 99.9%
<b>Ref No:GM10</b>	

CryoFit<sup>®</sup> is a trademark of Gatamed

**Gatamed Ltd.Sti.**

Asagi Eglence mah. Tevfik Saglam cad. No:7/3 06010, Kecioren, Ankara, TURKIYE  
Tel:+90 (532) 261 89 26 Fax:+90 (312) 321 96 98 sureyya@gatamed.com.tr www.gatamed.com